

new!

BAILEY

Bariatric Parallel Bars

Bariatric Parallel Bar Features:

- Adjustable height, 26" to 44", adjustable width, +/- 6".
- 10' Long, 1.5" diameter handrails.
- Offset handrail support provides unobstructed hand clearance.
- Durable chrome plated, tubular steel construction.
- ***Six post design allows 500 lb patient capacity!***
- **Floor mounting and assembly required.**
- Ship weight 100 lbs.

NOTE:

Floor mounting hardware is provided by others
to suit installation requirements.



**Bariatric Adjustable Height
and Width Parallel Bars**
[10' Long, Floor Mounted]

See our full catalog of products at
www.baileymfg.com