

BAILEY

MAT TABLES

STATIONARY & WALL MOUNTED MAT TABLES

UPHOLSTERED

Mat tables have been specially designed to provide an ample exercise area to the height most convenient to wheelchair and ambulatory people. Upholstered top table may be joined to form a larger exercise area.

MODEL 454 (4' X 7' X 20")



MODEL 456 & 457 (5' X 7' X 20" & 6' X 8' X 20")



MODEL 4520 BARIATRIC (4' X 7' X 28")



- 6 leg and special frame design allows 1000 pound patient capacity
- 2" upholstered top - US 387 Camel (standard) - other Bailey upholstery colors available
- Walnut finish available


RAISED RIM

Raised Rim Mat Tables feature a 1 1/4" rounded raised rim. The table surface is resin free to prevent damage to the accessory mats. Mats for Raised Rim Mat Tables are available as an option.

MODEL 450 (4' X 7' X 18") accessory 2" mat - **MODEL 51**

MODEL 452 (5' X 7' X 18") accessory 2" mat - **MODEL 53**

MODEL 453 (6' X 8' X 18") accessory 2" mat - **MODEL 54**



WALL MOUNTED- PLATFORM

This efficient, space saving, upholstered top mat platform locks into position against the wall when not in use. Extending out a mere 7" when folded, it leaves the surrounding floor space clear for other use.

MODEL 458 (4' X 7' X 20")

- 2" upholstered top - Herculite (R) #20 Navy (standard)
- mounting hardware not included
- extends 7" from wall when folded

ACCESSORIES


Made to fit Raised Rim Mat Tables, the extra big sizes are ideal for floor exercises. Two handles on one long side make lifting and hanging very easy. Zipper closure.

MODEL 50 (3' X 7' X 2")

MODEL 51 (4' X 7' X 2")


MODEL 53 CENTER SEAM (5' X 7' X 2")

MODEL 54 CENTER SEAM (6' X 8' X 2")



MODEL 741 MAT HANGER

- hang up to four 2" thick mats, used in pairs
- when not in use, the hook swings to the side
- brackets and hooks are made of strong electro-plated steel and are mounted on wooden bases
- mounting hardware not included



MODEL 742 MOBILE MAT TRUCK

- storage of up to eight 2" thick mats 6' X 8' or smaller

BAILEY MANUFACTURING CO.

P.O. BOX 130 LODI, OHIO 44254-1030

IN USA: PH 800-321-8372 FAX: 800-224-5390

OUTSIDE USA: PH 330-948-1080 FAX 330-948-4439

WWW.BAILEYMFG.COM



bailey_mfg_lodi



Bailey Manufacturing Company

CONTACT: